



PRODUCT CATALOGUE



## Your partner in semen analysis and Assisted Reproductive Technologies

Founded in 1992, FertiPro is a reliable, independent manufacturer of in vitro diagnostics and cell culture media used in assisted reproductive technologies and the diagnosis of male infertility.

### Product Range

We manufacture in vitro diagnostic kits for the determination of neutral alpha-glucosidase, fructose, anti-sperm antibodies, sperm vitality, white blood cells and sperm morphology in human semen samples. Our media range includes cell culture media for use in all stages of IVF and ICSI procedures as well as pre-washed, paraffin oils, a ready-to-use hyaluronidase solution, PVP solution, density gradients for sperm selection and media for cryopreservation of human oocytes, spermatozoa and embryos.

### Manufacturing

Manufacturing takes place in our own, state-of-the-art, fully environmentally controlled, production facility in Beernem, Belgium.

### Certification & quality assurance


FertiPro is an ISO 13485 certified company, that participates in the MDSAP program. We have product registrations in Europe, Brazil, Canada, United States of America and in increasingly more countries worldwide. A complete list of all distributors can be found on our website. All of our diagnostic kits and ART media are CE marked either as in vitro diagnostics or as medical devices. Certificates of analysis for all our products are available upon request.

### Contract manufacturing


Contact us for more information regarding contract manufacturing for your company.



## Semen analysis: product range diagnostics



## Assisted Reproductive Technologies: ART media





## Semen analysis: product range diagnostics



---

### Morphology

08 Spermac Stain

---

### Immunology

09 SpermMar Test IgA & IgG

---

### Vitality

10 VitalScreen

11 Hypo-osmotic Swelling Test

---

### Presence of white blood cells

12 LeucoScreen Plus

---

### Function of accessory glands

13 Fructose Test

14 EpiScreen Plus

# Spermac Stain

## Sperm morphology stain

### Regulatory

Europe: CE-marked (IVDR)  
 USA: registered – Canada: registered  
 Brazil: registered – Australia: registered  
 Other regions: information available upon request

### Storage

+25°C  
 Store between 2-25°C

Stable until expiry date specified on the label (approximately 27 months)

### Product order codes

SPS050 Spermac Stain – 4x 50ml  
 SPS250 Spermac Stain – 4x 250ml

### Product instruction video



An instruction film on the use of Spermac Stain is available on YouTube.




### General information

Sperm morphology analysis is one of the basic semen examinations performed in the diagnosis and management of male infertility. Spermac Stain is an in vitro diagnostic staining kit consisting of a fixative and 3 staining solutions for human spermatozoa. Staining facilitates distinction between morphologically normal and abnormal spermatozoa, and enhances visualization of different parts of the sperm cell (head, acrosome, equatorial region, midpiece, tail).

For professional use only

### Intended use

Spermac Stain is a qualitative, non-automated, diagnostic kit for professional use for the staining of human spermatozoa.



### General information

As sperm does not come into contact with the blood circulation, the male reproductive system contains no antisperm antibodies in normal conditions. However, when the blood-testis barrier is breached, the immune system can detect mature sperm as antigenic and form antisperm antibodies that cause sub- or infertility. Anti-sperm antibodies belong to two immunological classes: immunoglobulin (Ig)A and IgG antibodies, and can be present in the semen sample as well as in male blood serum. Antisperm IgA and IgG antibodies are clinically associated with immunological infertility, and screening can therefore provide help in assessing the male fertility.

For professional use only

### Intended use

The SpermMar Test IgA and IgG are both semi-quantitative, non-automated diagnostic kits for detecting antisperm antibodies of respectively the IgA or IgG class on spermatozoa in human semen (or serum when using SpermMar Test IgG).

# SpermMar Test IgA & IgG

## Test for the detection of anti-sperm antibodies

### Regulatory

Europe: CE-marked (IVDD) – USA: registered  
 Canada: Health Canada License  
 Brazil: registered – Australia: registered  
 Other regions: information available upon request

### Storage

+2°C Store between 2-8°C

+8°C

SpermMar Test IgA: 12 months from date of manufacture  
 SpermMar Test IgG: 18 months from date of manufacture

### Product order codes

SPMA_S	SpermMar Test IgA single kit	50 tests
SPMA_C	SpermMar Test IgA complete kit*	50 tests
SPMG_S	SpermMar Test IgG single kit	50 tests
SPMG_C	SpermMar Test IgG complete kit*	50 tests
SPMG_P	SpermMar Test IgG Positive Control	50 tests
SPMG_N	SpermMar Test IgG Negative Control	50 tests

### Product instruction video



Instruction films on the use of SpermMar Test IgA & IgG are available on YouTube.

\* The complete kit contains pipettes and glassware to perform the test.

# VitalScreen

Sperm vitality test

## Regulatory

Europe: CE-marked (IVDR) – USA: registered – Canada: registered  
Brazil: registered – Australia: registered  
Other regions: information available upon request

## Storage

+25°C  
+2°C Store between 2-25°C  
24 months from date of manufacture

## Product order codes

VITAL VitalScreen 200 tests

## Product instruction video



An instruction film on the use of VitalScreen is available on YouTube.



## General information

VitalScreen uses eosin-nigrosin staining to assess the percentage of live spermatozoa. The technique is based on the principle that dead cells take up eosin through their cell membrane, and as a result stain red. The nigrosin provides a dark background which makes it easier to evaluate the slides. Sperm vitality should be determined when less than 40% of spermatozoa are motile. In samples with poor motility, it is important to discriminate between immotile dead sperm and immotile live sperm.

*For professional use only*

## Intended use

VitalScreen is a semi-quantitative, non-automated diagnostic test to evaluate the vitality of spermatozoa in a semen sample by using an eosin-nigrosin staining.

## General information

The Hypo-osmotic Swelling Test (HOS Test) is an in vitro diagnostic for professional use to evaluate the vitality of the spermatozoa in a semen sample. In contrast to dead spermatozoa, living cells have intact cell membranes which allow regulated water transport in hypo-osmotic conditions which will result in swelling or curling of the sperm tail upon incubation in HOST medium.

*For professional use only*



## Intended use

The HOS Test is a semi-quantitative, non-automated diagnostic test to evaluate the vitality (membrane function) of spermatozoa in a semen sample. The HOS Test should not be used for the selection of sperm in ART procedures such as intra-cytoplasmic sperm injection (ICSI).

# Hypo-osmotic Swelling Test

Sperm vitality test

## Regulatory

Europe: CE-marked (IVDD) – Canada: Health Canada License  
Brazil: registered – Australia: registered – USA: registered  
Other regions: information available upon request

## Storage

+8°C  
+2°C Store between 2-8°C  
12 months from date of manufacture

## Product order codes

HOST Hypo-osmotic Swelling Test 5x 20ml

## Product instruction video



An instruction film on the use of the Hypo-osmotic Swelling Test is available on YouTube.

# LeucoScreen Plus

Improved leucocyte staining method

## Regulatory

Europe: CE-marked (IVDD) – USA: registered  
Australia: registered – Brazil: registered – Canada: registered  
Other regions: information available upon request

## Storage

+25°C  
+2°C  
Store between 2-25°C  
12 months from date of manufacture

## Product order codes

LEUCO\_PLUS LeucoScreen Plus 40 days of analysis

## Product instruction video



An instruction film on the use of LeucoScreen Plus is available on YouTube. Training form can be downloaded from [www.fertipro.com](http://www.fertipro.com)



## General information

Most human ejaculates contain leucocytes and the predominant form of leucocytes in human semen are peroxidase-positive granulocytes. Excessive presence of these cells (leucocytospermia) may indicate the existence of reproductive tract infection. Leucocytospermia may also be associated with defects in the semen profile (reduction in sperm motility and DNA integrity, raise of sperm viscosity as well as loss of sperm function as a result of oxidative stress, and/or secretion of cytotoxic cytokines by these white blood cells). Although leucocytospermia is not an absolute indication of infertility, this condition is observed on average in 10 to 20% of all infertile men.

*For professional use only*

## Intended use

LeucoScreen Plus is a semi-quantitative, non-automated, histochemical and diagnostic kit for the determination of peroxidase-positive white blood cells in human semen.



## General information

The secretion of the seminal vesicles constitutes the main fraction of the ejaculate. Seminal fructose concentration can be used as a marker of the seminal vesicular function. Low fructose in semen is characteristic for partial or complete ejaculatory duct obstruction, bilateral congenital absence of the vas deference, partial retrograde ejaculation and androgen deficiency.

*For professional use only*

## Intended use

The Fructose Test is a quantitative, non-automated, photometric and diagnostic kit for the measurement of fructose in human semen or seminal plasma (fresh or frozen).

# Fructose Test

Kit for the determination of fructose in seminal plasma

## Regulatory

Europe: CE-marked (IVDD) – Australia: registered  
Other regions: information available upon request

## Storage

+8°C  
+2°C  
Store between 2-8°C  
12 months from date of manufacture

## Product order codes

FRUCTOSE Fructose Test 96 tests

## Product instruction video



An instruction film on the use of the Fructose Test is available on YouTube

# EpiScreen Plus

Kit for the determination of neutral alpha-glucosidase in seminal plasma

## Regulatory

Europe: CE-marked (IVDD) – Australia: registered  
Other regions: information available upon request

## Storage

+2°C Store between +2°C and +8°C  
24 months from date of manufacture

## Product order codes

EPI\_PLUS EpiScreen Plus – 25 tests

## Product instruction video



An instruction film on the use of EpiScreen Plus is available on YouTube.



## General information

EpiScreen Plus can be used to determine the neutral alpha-glucosidase activity in semen (plasma), an enzyme which is mainly secreted by the epididymis. The activity of this enzyme is a reliable marker for epididymis function in patients with (very) low sperm concentration or azoospermic patients, having a normal androgen blood level:

- very low activity indicates a bilateral obstruction between the epididymis and the ejaculatory duct.
- low activity may reflect partial obstruction of the epididymis.
- normal enzyme activity is expected when there is an obstruction above the area in which the enzyme is secreted or in cases of non-obstructive azoospermia (testicular dysfunction).

For professional use only

## Intended use

EpiScreen Plus is a quantitative, non-automated, photometric and diagnostic kit for detecting neutral alpha-glucosidase in human semen or seminal plasma.





# Assisted Reproductive Technologies: ART media



## Sperm processing

### Washing/swim-up

- 18 FertiCult Flushing medium
- 20 Sil-Select Plus Sperm Washing/ Insemination medium
- 28 GAIN medium

### Density gradient separation

- 20 Sil-Select Plus
- 22 Sil-Select STOCK

### Sperm processing for ICSI

- 23 10% PVP in FertiCult Flushing medium



## Intrauterine insemination

- 18 FertiCult Flushing medium
- 20 Sil-Select Plus Sperm Washing/ Insemination medium
- 28 GAIN medium



## Oocyte handling

### Washing

- 18 FertiCult Flushing medium
- 24 FertiCult Aspiration medium

### Pick-up

- 24 FertiCult Aspiration medium

### Denudation

- 25 Hyaluronidase in FertiCult Flushing medium

### Oocyte incubation

- 26 FertiCult IVF medium
- 28 GAIN medium



## Fertilization

- 26 FertiCult IVF medium
- 28 GAIN medium
- 30 FertiCult Mineral Oil
- 31 FertiCult High Viscosity Oil



## Embryo handling

### Embryo washing

- 18 FertiCult Flushing medium

### Embryo culture

- 28 GAIN medium

### Culture overlay

- 29 FertiCult Mineral Oil
- 29 FertiCult High Viscosity Oil



## Embryo transfer

- 18 FertiCult Flushing medium
- 28 GAIN medium



## Cryopreservation media

### Sperm cryopreservation

- 32 SpermFreeze
- 32 SpermFreeze SSP
- 33 SpermFreeze Box

### Oocyte cryopreservation

- 34 FertiVit Cooling / FertiVit Warming

### Embryo cryopreservation

- 35 VitriFreeze / VitriThaw
- 36 FertiVit Cooling / FertiVit Warming
- 37 VitriFreeze ES / VitriThaw ES

# FertiCult Flushing medium

HEPES buffered cell culture medium for use in human ART



## Intended use

FertiCult Flushing medium is intended for in vitro procedures involving human gametes (sperm and oocytes), including washing of gametes, sperm swim-up procedures, intra-uterine insemination (IUI) of the spermatozoa and sperm injection during intracytoplasmic sperm injection (ICSI). FertiCult Flushing medium can also be used for human embryo washing and holding, and for embryo transfer (ET).

*For professional use only*



## Product specifications

- Chemical composition
- pH: 7.30-7.90 (release criteria: 7.30-7.60)
- Osmolality: 270-290 mOsm/kg
- Endotoxin test by Limulus Amebocyte Lysate (LAL) methodology (USP <85>): < 0.25 EU/ml
- Sterility test by the current Ph. Eur. 2.6.1 / USP <71>: No growth
- One-cell mouse embryo assay (% blastocysts after 96 hours, after 60 minutes exposure to test medium): ≥ 80%
- Human sperm survival assay (% motility compared with control after 24 hours exposure to test medium): ≥ 80%
- Use of Ph Eur or USP grade products if applicable

## Regulatory

Europe: CE marked (MDR, Notified Body number 2797)

Canada: Health Canada Licence – Brazil: registered

Other regions: information available upon request

## Storage

+2°C  Store between 2-25°C

+2°C  Store between 2-8°C (for the gentamicin containing products)

 18 months from date of manufacture

**STERILE A** Sterilized using aseptic processing techniques

 Keep away from (sun)light

## Product order codes

FLUSH020	FertiCult™ Flushing medium	5x 20ml
FLUSH050	FertiCult™ Flushing medium	5x 50ml
FLUSH100	FertiCult™ Flushing medium	3x 100ml
FLUSH020PHR	FertiCult™ Flushing medium with phenol red	5x 20ml
FLUSH050PHR	FertiCult™ Flushing medium with phenol red	5x 50ml
FLUSH100PHR	FertiCult™ Flushing medium with phenol red	3x 100ml
FLUSH500PHR	FertiCult™ Flushing medium with phenol red	1x 500ml
FLUSH020PRG	FertiCult™ Flushing medium with phenol red and gentamicin	5x 20ml
FLUSH050PRG	FertiCult™ Flushing medium with phenol red and gentamicin	5x 50ml
FLUSH100PRG_G	FertiCult™ Flushing medium with phenol red and gentamicin	3x 100ml
FLUSH500PRG_G	FertiCult™ Flushing medium with phenol red and gentamicin	1x 500ml
WASH005	Sil-Select Plus™ Sperm Washing/ Insemination medium	25x 5 ml
WASH020	Sil-Select Plus™ Sperm Washing/ Insemination medium	5x 20 ml
WASH100G	Sil-Select Plus™ Sperm Washing/ Insemination medium with gentamicin	1x 100ml

# Sil-Select Plus

## Sperm preparation media



### Intended use

Sil-Select Plus is a ready-to-use gradient system for semen preparation.

The following gradients are available:

- Sil-Select Plus Upper layer (45%)
- Sil-Select Plus 80%
- Sil-Select Plus Lower layer (90%)

The Sil-Select Plus system can be used for semen preparation for Intra Uterine Insemination (IUI), In Vitro Fertilization (IVF) and IntraCytoplasmic Sperm Injection (ICSI).

Sil-Select Plus Sperm Washing/Insemination medium is intended for use in Assisted Reproductive Technologies, more specifically for sperm washing/holding and for the production of density gradients. The medium can be used for IUI.

*For professional use only*



### Product instruction video



An instruction film on the use of Sil-Select Plus is available on our website.

### Product specifications

#### Sil-Select Plus Upper layer/Lower layer/80% media:

- Chemical composition
- pH criteria: 7.20-7.90 (release criteria: 7.20-7.60)
- Osmolality: 280-310 mOsm/kg (Upper layer)  
290-330 mOsm/kg (Lower layer)  
290-330 mOsm/kg (80%)
- Density: 1.0500-1.0700 (Upper layer)  
1.1050-1.1150 g/ml (Lower layer)  
1.097-1.107 g/ml (80%)
- Endotoxin test by Limulus Amebocyte Lysate (LAL) methodology (USP <85>): < 0.5EU/ml
- Sterility test according to the current Ph. Eur. 2.6.1/USP <71>: No growth
- Human sperm survival assay (% motility compared with control after 4 hours exposure to test medium): ≥ 80%
- Human sperm survival assay (% motility compared with control after 24 hours exposure to test medium): ≥ 75%
- Not MEA tested
- Use of Ph Eur or USP grade products if applicable

#### Sil-Select Plus Sperm Washing/Insemination medium:

- Chemical composition
- pH criteria: 7.30-7.90 (release criteria: 7.30-7.60,)
- Osmolality: 270-290 mOsm/kg
- Sterility test according to the current Ph. Eur. 2.6.1/USP <71>: No growth
- Endotoxin test by Limulus Amebocyte Lysate (LAL) methodology (USP <85>): < 0.25EU/ml
- One-cell Mouse Embryo Assay (blastocysts after 96 hours, after 60 minutes exposure to test medium): ≥ 80%
- Human Sperm Survival Assay (% motility compared with control after 24 hours exposure to test medium): ≥ 80%
- Use of Ph Eur or USP grade products if applicable

### Regulatory

Europe: CE-marked (MDR, Notified Body number 2797)

Canada: Health Canada License – Brazil: Registered

USA: US FDA-cleared – Other regions: Information available upon request

### Product order codes

SIP008	Sil-Select Plus™ Lower layer	8x 2.5ml
	Sil-Select Plus™ Upper layer	8x 2.5ml
	Sil-Select Plus™ Sperm Washing/Insemination medium	8x 5ml
SIP016	Sil-Select Plus™ Lower layer	2x 20ml
	Sil-Select Plus™ Upper layer	2x 20ml
	Sil-Select Plus™ Sperm Washing/Insemination medium	5x 20ml
SIP020LO	Sil-Select Plus™ Lower layer	1x 20ml
SIP020UP	Sil-Select Plus™ Upper layer	1x 20ml
SIP050	Sil-Select Plus™ Lower layer	1x 50ml
	Sil-Select Plus™ Upper layer	1x 50ml
SIP100	Sil-Select Plus™ Lower layer	1x 100ml
	Sil-Select Plus™ Upper layer	1x 100ml
SIP80_100	Sil-Select Plus™ 80%	1x 100ml

*Single bottles can be ordered separately*

### Storage



Store between  
2-8°C



18 months from date  
of manufacture



Keep away  
from (sun)light



Sterilized using aseptic  
processing techniques

# Sil-Select STOCK

Sperm preparation media



# 10% PVP in FertiCult Flushing medium

Viscous medium for use in ICSI procedures

## Product specifications

- Chemical composition
- pH criteria: 7.20-7.90 (Release criteria: 7.20-7.60)
- Osmolality: 300-330 mOsm/kg
- Density: 1.1150-1.1250 g/ml
- Endotoxin test by Limulus Amebocyte Lysate (LAL) methodology (USP <85>): < 0.5EU/ml
- Sterility test according to the current Ph. Eur. 2.6.1/USP <71>: No growth
- Human sperm survival assay (% motility compared with control after 4 hours exposure to test medium): ≥ 80%
- Human sperm survival assay (% motility compared with control after 24 hours exposure to test medium): ≥ 75%
- Not MEA tested
- Use of Ph Eur or USP grade products if applicable

## Regulatory

Europe: CE-marked (MDR, Notified Body number 2797)  
 USA: US FDA Cleared  
 Canada: Health Canada License  
 Brazil: Registered  
 Other regions: Information available upon request

## Storage



Store between 2-8°C



18 months from date of manufacture



Sterilized using aseptic processing techniques



Keep away from (sun)light



## Intended use

Sil-Select Stock is a stock solution for semen preparation. It is an isotonic gradient for semen preparation with a density of approximately 1.12 g/ml. Sil-Select Stock can be used for semen preparation for Intra Uterine Insemination (IUI), In Vitro Fertilization (IVF) and IntraCytoplasmatic Sperm Injection (ICSI).

For professional use only



## Intended use

10% PVP in FertiCult Flushing medium is a viscous medium containing polyvinylpyrrolidone (PVP) used for sperm preparation for Intracytoplasmic Sperm Injection (ICSI). These procedures require the capture of individual sperm cells in a glass pipette for injection into the oocyte and this is facilitated by first immobilizing the sperm by placing them in a viscous medium like 10% PVP in FertiCult Flushing medium prior to nicking the tail to immobilize sperm completely.

For professional use only

## Product specifications

- Chemical composition
- pH: 7.20-7.90 (release criteria: 7.20-7.60)
- Osmolality: 280-340 mOsm/kg (release criteria: 280-320 mOsm/kg)
- Viscosity: >170 cP at 25 °C
- Endotoxin test by Limulus Amebocyte Lysate (LAL) methodology (USP <85>): < 0.5EU/ml
- Sterility test by the current Ph. Eur. 2.6.1 / USP <71>: No growth
- One-cell mouse embryo assay (% blastocysts after 96 hours, after 10 minutes exposure time to test medium): ≥ 80%
- Human sperm survival assay (% motility compared with control, after 60 minutes exposure time to test medium): ≥ 80%
- Use of Ph Eur or USP grade products if applicable

## Regulatory

Europe: CE-marked (MDR, Notified Body number 2797)  
 Canada: Health Canada License  
 Brazil: Registered  
 Other regions: Information available upon request

## Storage



Store between 2-8°C



9 months from date of manufacture



Sterilized using aseptic processing techniques



Keep away from (sun)light

## Product order codes

PVP0.2S	10% PVP in FertiCult Flushing medium (screw cap)	5x 0.2ml
PVP1	10% PVP in FertiCult Flushing medium	1x 1ml

## FertiCult Aspiration medium

HEPES buffered cell culture medium for human oocyte pick-up and follicular flushing



## Hyaluronidase in FertiCult Flushing medium

Hyaluronidase solution

### Product specifications

pH (ambient air)	7.30-7.90 (release criteria: 7.30-7.60)
Osmolality	270-290 mOsm/kg
Sterility	No growth
Endotoxin	< 0.25 EU/ml
Mouse embryo assay*	≥ 80% blastocysts formed after 96h incubation
Use of Ph Eur or USP grade products if applicable	

### Regulatory

Europe: CE-marked (Notified Body 0344)  
Brazil: registered  
Other regions: information available upon request

### Product order codes

ASPI020	FertiCult Aspiration medium	5x 20ml
ASPI100	FertiCult Aspiration medium	5x 100ml

### Storage



Store between  
2-25°C



18 months from date  
of manufacture



Keep away  
from sunlight



### Intended use

FertiCult Aspiration medium is a ready to use formulation for oocyte pick-up and follicular flushing. FertiCult Aspiration medium contains heparin (pharmaceutical grade) and HEPES to stabilize the pH under air, no CO<sub>2</sub> incubation is required. Product is sterilized using aseptic processing techniques.

For professional use only



### Intended use

Hyaluronidase in FertiCult Flushing medium contains derived 80 IU/ml hyaluronidase in FertiCult Flushing medium. The solution is used in the oocyte denudation process. Hyaluronidase digests the hyaluronic acid between the cumulus cells, which makes it easier to remove the cumulus mechanically. Product is sterilized using aseptic processing techniques.

For professional use only

### Product specifications

pH (ambient air)	7.30-7.90 (release criteria: 7.30-7.60)
Osmolality	270-290 mOsm/kg
Sterility	No growth
Endotoxin	< 1.00 EU/ml
Mouse embryo assay*	≥ 80% blastocysts formed after 96h incubation
Use of Ph Eur or USP grade products if applicable	

### Regulatory

Europe: CE-marked (Notified Body 0344)  
Brazil: registered  
Other regions: information available upon request

### Product order codes

HYA001	Hyaluronidase in FertiCult Flushing medium	5x 1ml
HYA010	Hyaluronidase in FertiCult Flushing medium	1x 10ml

### Storage



Store between  
2-8°C



12 months from date  
of manufacture



Keep away  
from sunlight

# FertiCult IVF medium

## Fertilization medium



### Intended use

FertiCult IVF medium is intended for washing and holding of human oocytes, performing fertilization by In Vitro Fertilization (IVF) or IntraCytoplasmic Sperm Injection (ICSI) (until 2PN).

*For professional use only*



### Product specifications

- Chemical composition
- pH: 7.20-7.50 (37 °C - 6% CO<sub>2</sub>)
- Osmolality: 270-290 mOsm/kg
- Sterility test by the current Ph. Eur. 2.6.1. / USP <71>: No growth
- Endotoxin test by Limulus Amebocyte Lysate (LAL) methodology (USP <85>): < 0.25EU/ml
- Mouse Embryo Assay (% blastocysts after 96h after IVF in test medium): ≥ 80%
- Use of Ph Eur or USP grade products if applicable

### Product order codes

FECU020	FertiCult™ IVF medium	5x 20ml
FECU050	FertiCult™ IVF medium	5x 50ml
FECU100	FertiCult™ IVF medium	3x 100ml
FECU020PHR	FertiCult™ IVF medium with phenol red	5x 20ml
FECU050PHR	FertiCult™ IVF medium with phenol red	5x 50ml
FECU100PHR	FertiCult™ IVF medium with phenol red	3x 100ml
FECU020PHR_G	FertiCult™ IVF medium with phenol red and gentamicin	5x 20ml
FECU050PHR_G	FertiCult™ IVF medium with phenol red and gentamicin	5x 50ml
FECU100PHR_G	FertiCult™ IVF medium with phenol red and gentamicin	3x 100ml

### Regulatory

Europe: CE-marked MDR (Notified Body number 2797)

Brazil: registered

Other regions: information available upon request

### Storage



Store between  
2-25°C



12 months  
from date of manufacture



Sterilized using aseptic  
processing techniques



Store between 2-8°C  
(for the gentamicin  
containing products)



Keep away  
from (sun)light

# GAIN medium

Cell culture media for use in human ART



## Product specifications

- Chemical composition
- pH: 7.20-7.45 (37 °C, 5% CO<sub>2</sub>)
- Osmolality: 270-290 mOsm/kg
- Sterility test by the current Ph. Eur. 2.6.1/ USP <71>: No growth
- Endotoxin test by Limulus Amebocyte Lysate (LAL) methodology (USP <85>): < 0.25EU/ml
- One-cell mouse embryo assay (% blastocysts after 96 hours culture in test medium): ≥ 80%
- Human sperm survival assay (% motility compared with control, after 24 hours exposure to test medium at room temperature, no CO<sub>2</sub> incubation): ≥ 80%
- Use of Ph Eur or USP grade products if applicable

## Regulatory

Europe: CE-marked (MDR, Notified Body Number 2797)  
 Brazil: Registered  
 Other regions: Information available upon request

## Storage



STERILE A

Sterilized using aseptic processing techniques



Keep away from sunlight

## Product order codes

GAIN010	GAIN™ medium	5 x 10ml
GAIN020	GAIN™ medium	5 x 20ml



## Intended use

GAIN medium is a single-step cell culture medium for use with human embryos and gametes. GAIN medium can additionally be used in all the procedures indicated below:

- For semen washing and Intra Uterine Insemination (IUI).
- For oocytes handling/incubation in preparation of, or during fertilization by In Vitro Fertilization (IVF)/ Intra Cytoplasmic Sperm Injection (ICSI).
- For embryo culture from day 1 to expanded blastocyst stage.
- For embryo transfer (ET).

For professional use only

**NEW**

# FertiCult High Viscosity Oil



Experience the comfort of using FertiCult High Viscosity Oil with significant higher viscosity compared to FertiCult Mineral Oil!

- ✓ **Make it simple. Preparing and handling dishes will never have been that easy!**
- ✓ **Feel more relaxed. Sometimes, manipulations outside the incubator need more time. The viscous oil layer guarantees a perfectly stable pH for at least 30 minutes!**

FertiCult Mineral Oil and FertiCult High Viscosity Oil are produced and tested according to the highest quality standards and CE marked according to the EU-MDR.

# FertiCult Mineral Oil - FertiCult High Viscosity Oil

Pre-washed mineral oils



**NEW**

## Intended use

FertiCult Mineral Oil and FertiCult High Viscosity Oil are ready to use, pre-washed mineral oils for use in In Vitro Fertilization (IVF), Intracytoplasmic Sperm Injection (ICSI) and related Assisted Reproductive Technologies (ART).

FertiCult High Viscosity Oil is more viscous and has higher density values compared to FertiCult Mineral Oil. The difference in viscosity (and density) is reflected in the practical use of the product and is dependent on the procedures used by a certain laboratory. In this way, laboratories can choose between both types of oil in order to guarantee the most optimal culture system according to their procedures.

For professional use only

## Product specifications

	FertiCult Mineral Oil	FertiCult High Viscosity Oil
Chemical composition	✓	✓
Density:	0.83-0.86g/ml	0.86-0.88g/ml
Viscosity:	< 60cP at 30 °C	95-140cP at 30 °C
Endotoxin test by Limulus Amebocyte Lysate (LAL) methodology (USP <85>):	< 0.1 EU/ml	< 0.1 EU/ml
Peroxide Value (POV):	<0.1 mEq/kg	< 0.1 mEq/kg
Sterility test by the current Ph. Eur. 2.6.1/ USP <71>:	No growth	No growth
One-cell Mouse embryo assay %:	≥ 80% (blastocysts after 120 hours in standard mouse IVF medium covered with test item)	≥ 80% (blastocysts after 120 hours in standard mouse IVF medium covered with test item)
Use of pharmaceutical grade mineral (paraffin)oil	✓	✓

## Regulatory

Europe: CE-marked (MDR, Notified Body Number 2797)  
 USA: US FDA Cleared  
 Brazil: registered (FertiCult Mineral Oil)  
 Other regions: information available upon request

## Product order codes

(glass bottle)  
 MINOIL050 FertiCult™ Mineral Oil 1 x 50ml  
 MINOIL100 FertiCult™ Mineral Oil 1 x 100ml  
 MINOIL500 FertiCult™ Mineral Oil 1 x 500ml  
 (PETG bottle)  
 MINOIL050P FertiCult™ Mineral Oil 1 x 50ml  
 MINOIL100P FertiCult™ Mineral Oil 1 x 100ml  
 MINOIL500P FertiCult™ Mineral Oil 1 x 500ml

## Product order codes

(glass bottle)  
 HVOIL050 FertiCult™ High Viscosity Oil 1 x 50ml  
 HVOIL100 FertiCult™ High Viscosity Oil 1 x 100ml  
 HVOIL500 FertiCult™ High Viscosity Oil 1 x 500ml

## Storage





# SpermFreeze / SpermFreeze SSP

Medium for freezing human spermatozoa



## Product specifications

- Chemical composition
- pH 7.20 – 7.90 (Release criteria: 7.20 – 7.60)
- Sterility test by the current Ph. Eur. 2.6.1 / USP <71>: No growth
- Endotoxins (USP <85>): < 0.25EU/ml
- Human sperm survival assay (% motility compared with control after 4 hours exposure to test medium): ≥ 80%
- Not MEA tested
- Use of Ph Eur or USP grade products if applicable

## Regulatory

Europe: CE-marked (MDR, Notified Body number 2797)  
USA: US FDA Cleared  
Canada: Health Canada license (only SpermFreeze)  
Brazil: Registered  
Other regions: Information available upon request

## Product order codes

### SpermFreeze

SPF05 SpermFreeze™ 25 x 5ml  
SPF SpermFreeze™ 5 x 20ml

### SpermFreeze SSP

SSP001 SpermFreeze™ SSP 5 x 1ml

## Storage




## Intended use

SpermFreeze and SpermFreeze SSP are media for cryopreservation of human sperm for further use in Assisted Reproductive Technologies.

For professional use only

**NEW**




# SpermFreeze Box

- ✓ Freezing tool to standardize your sperm freezing protocols
- ✓ Ideal height above liquid nitrogen level guaranteed thanks to the SpermFreeze Floater
- ✓ Robust Freezing Box with a thick wall to work safe with liquid nitrogen
- ✓ Manufactured from EPP (Expanded Poly Propylene)
- ✓ Multiple use, washable (dish washer) and sterilisable

## Protocol

- 1 Fill the box with 1 liter of liquid nitrogen
- 2 Place your sealed straws on the floater
- 3 Deposit the floater on the liquid nitrogen
- 4 Close the box with the lid
- 5 Wait for at least 15 minutes
- 6 Plunge the straws in the liquid nitrogen
- 7 Store them in your nitrogen storage container for long term preservation

Product order code: **SPF\_BOX** SpermFreeze Box (342 x 264 x 188mm)



# VitriFreeze – VitriThaw

Kits for vitrifying and warming of human blastocysts



## Intended use

VitriFreeze and VitriThaw are a set of media for vitrification and thawing of human embryos (morula till expanded blastocyst stage).

*For professional use only*

## Product specifications

- Chemical composition
- pH: 7.20-7.40
- Osmolality (mOsm/kg):
  - VitriFreeze Pre-incubation: 270-295 (release criteria: 270-290)
  - VitriThaw Thawing 1: 805-865 (release criteria: 805-850)
  - VitriThaw Thawing 2: 535-565
  - VitriThaw Thawing 3: 405-435
  - VitriThaw Thawing 4: 270-295 (release criteria: 270-290)
- Sterility test by the current Ph. Eur. 2.6.1. / USP <71>: no growth
- Endotoxins test by Limulus Lysate methodology (USP <85>): < 0.25 EU/ml
- One-Cell Mouse embryo assay (% blastocysts stage after 96 hours following a sequential exposure to the cooling and warming media) ≥ 80%
- Use of Ph Eur or USP grade products if applicable

## Product order codes

VF\_KIT1 VitriFreeze™ kit  
 VT\_KIT1 VitriThaw™ kit

## Storage

+2°C to +8°C  
 Store between 2-8°C

12 months from date of manufacture

## Regulatory

Europe: CE-marked (Notified Body 0344)  
 USA: US FDA Cleared – Other regions:  
*information available upon request*

STERILE A

Sterilized using aseptic processing techniques

Keep away from sunlight

## FertiVit Cooling – FertiVit Warming

Kits for vitrification and warming of human oocytes and embryos up till blastocyst stage



### Product specifications

- Chemical composition.
- pH: 7,20-7,50 (release criteria: 7,20-7,40).
- Osmolality (mOsm/kg):
  - Pre-incubation / Warming 6: 270-295 (release criteria: 270-290)
  - Warming 3: 805-865 (release criteria: 805-850)
  - Warming 4: 535-565
  - Warming 5: 405-435
- Sterility test by the current Ph. Eur. 2.6.1. / USP <71>; No growth
- Endotoxins test by Limulus Amebocyte Lysate (LAL) methodology (USP <85>): < 0,25 EU/ml.
- One-Cell Mouse Embryo Assay (% blastocyst stage after 96 hours following a sequential exposure to the cooling and warming media): ≥ 80%.
- Use of Ph Eur or USP grade products if applicable.

### Product order codes

FVC\_KIT FertiVit™ Cooling kit  
FVW\_KIT FertiVit™ Warming kit

### Regulatory

Europe: CE-marked (Notified Body 0344)  
Brazil: registered  
Other regions: information available upon request

### Storage



Store between  
+2°C and +8°C

STERILE A

Sterilized using aseptic  
processing techniques



12 months from date  
of manufacture



Keep away  
from (sun)light



### Intended use

FertiVit Cooling and FertiVit Warming are a set of media for vitrification and warming of human oocytes and embryos (zygote till blastocyst).

For professional use only



### Intended use

VitriFreeze ES and VitriThaw ES are a set of media for vitrification and thawing of human embryos (zygote till expanded blastocyst stage).

For professional use only

## VitriFreeze ES – VitriThaw ES

Kits for vitrifying and warming of human embryos

### Product specifications

- Chemical composition
- pH: 7,20-7,40
- Osmolality (mOsm/kg):
  - VitriFreeze ES Pre-incubation medium: 270-295 (release criteria: 270-290)
  - VitriThaw ES Thawing medium 2: 805-865 (release criteria: 805-850)
  - VitriThaw ES Thawing medium 3: 535-565
  - VitriThaw ES Thawing medium 4: 405-435
  - VitriThaw ES Thawing medium 5: 270-295 (release criteria: 270-290)
- Sterility test by the current Ph. Eur. 2.6.1. / USP <71>; no growth
- Endotoxins test by Limulus Lysate (LAL) methodology (USP <85>): < 0,25 EU/ml
- One-cell Mouse Embryo Assay (% blastocysts stage after 96 hours following a sequential exposure to the cooling and warming media) ≥ 80%
- Use of Ph Eur or USP grade products if applicable

### Regulatory

Europe: CE-marked (Notified Body 0344)  
Brazil: registered  
Other regions: information available upon request

### Product order codes

VF\_KIT1\_ES VitriFreeze ES™ kit  
VT\_KIT1\_ES VitriThaw ES™ kit

### Storage



Store between  
+2°C and +8°C

STERILE A

Sterilized using aseptic  
processing techniques



12 months from date  
of manufacture



Keep away  
from (sun)light





# FertiPro

-  ISO 13485 certified
-  MDSAP certified
-  EU-MDR approved medical devices
-  EU-IVDR conform diagnostics



# FertiPro

FertiPro NV  
Industriepark Noord 32  
B-8730 Beernem  
Belgium

tel. +32 (0)50 79 18 05  
fax +32 (0)50 79 17 99  
info@fertipro.com  
www.fertipro.com

EDITION DATE 06 2022



Subscribe  
to our  
newsletter