


# SpermFreeze Box

**NEW**



- ✓ Freezing tool to standardize your sperm freezing protocols
- ✓ Ideal height above liquid nitrogen level guaranteed thanks to the SpermFreeze Floater
- ✓ Robust Freezing Box with a thick wall to work safe with liquid nitrogen
- ✓ Manufactured from EPP (Expanded Poly Propylene)
- ✓ Multiple use, washable (dish washer) and resistant to broad range of disinfectants

## FertiCult Flushing medium

HEPES buffered cell culture medium for use in human ART



# FertiPro



**NEW**

## High Viscosity Oil

Experience the comfort of using FertiCult High Viscosity Oil with significant higher viscosity compared to FertiCult Mineral Oil!

- ✓ *Preparing & handling dishes have never been that easy!*
- ✓ *The viscous oil layer guarantees a perfectly stable pH for at least 30 minutes!*

FertiCult Mineral Oil and FertiCult High Viscosity oil are produced and tested according to the highest quality standards and CE marked according to the EU-MDR.

## Discover our other products

