

HALOSPERM® G2 PRO

Ultra-Fast DNA Fragmentation Assessment



NEW



- DNA Fragmentation SCD Method
- Quickest Product in the Market
- 10 Tests/Kit



Halosperm® G2 PRO is the latest evolution in sperm DNA fragmentation testing, delivering the fastest processing time on the market—just 10 minutes. This new system offers IVF laboratories Halosperm® G2 kit, this third-generation system offers IVF laboratories an ultra-rapid, user-friendly, and odor-reduced workflow that requires no fume hood.

The kit includes optimized reagents and ready-to-use staining solutions, minimizing variability and enabling clear visualization under bright-field microscopy.

Key Features:

- **Fastest processing time** – available results in just 10 minutes
- **Simplified**, hood-free protocol for maximum lab flexibility
- Bright-field compatible staining for consistent, **high-quality visualization**
- Based on **proven SCD technology** for accurate DNA fragmentation analysis
- **CE marked** and available as an **IVD-certified** version
- **Cost effective**

Halosperm® G2 PRO is the new gold standard for rapid, easy, and reliable sperm DNA fragmentation testing.

HALOSPERM® G2 PRO

Determinations	10
Sample kit determination	2
Validity (months)	12
Stability - after kit opening (months)	9
Fume hood	Recommended
Stains	Brightfield Stain
Processment	Drop bottles
Processing time	10 min